
NeuroHawkTM Thrombectomy Device

Instructions for Use

MicroPort NeuroTech (Shanghai) Co., Ltd.

Content

Warnings	3
Product Description	3
Product Specifications/Use Guidelines	4
Intended Purpose	4
Indication For Use	4
Intended User	4
Target Population	4
Clinical Benefit	4
Contraindications	4
Potential Complications	5
Precautions	5
Accessories Required for Product Use	5
Operation Procedures	6
Disposal	8
Sterilization	9
Storage Conditions	9
Transportation Conditions	9
Shelf Life	9
Basic UDI-DI number	9
Legal Statement	9
After-sales Service	9
Legal Manufacturer/Manufacturing Site	9
After-sales Service Company Information	9
European Authorized Representative	10
Symbol Glossary	11

NeuroHawk™ Thrombectomy Device

Instructions for Use

Warnings

- This product has been sterilized by ethylene oxide.
- This product is for single-use only. Do not re-sterilize or reuse.
- Do not use if the packaging is found to have been opened, damaged, or a leak is detected, or the validity period has expired prior to use.
- This Instructions for Use must be read carefully before use.
- This product should be used by qualified professionals who have received necessary interventional training (especially for intracranial thrombectomy) to complete the endovascular procedures.
- The product should not be repositioned in the vessel without complete retrieval into the microcatheter and should be reopened only after complete retrieval into the microcatheter.
- When using this product for thrombus removal, it is recommended, based on the product characteristics, that retrieval attempts for a single thrombectomy stent should be less than or equal to 5, and that the retrieval attempts in the same vessel should be less than or equal to 5 in order to reduce the risk of vascular injury.
- Allergic reactions may occur with the use of this product in people who are allergic to Nitinol.
- After the device has been used, the product and packaging should be disposed of in accordance with the regulations of hospitals, administrative departments and local governments.

Product Description

NeuroHawk™ thrombectomy device is composed of a stent system and an introducer sheath (Figure 1), wherein the stent system consists of a self-expanding stent and a delivery wire.

The NeuroHawk™ thrombectomy device is supplied sterile for single use only. It is used in acute ischemic strokes caused by large intracranial vessel occlusion. The product adopts the minimally invasive interventional treatment method. Under the supervision of DSA, interventional access is first established by using a puncture sheath, guidewire, catheter, and microcatheter. The thrombectomy stent is then delivered within the microcatheter to the target site by the delivery wire of the thrombectomy device. The self-expandable thrombectomy stent is deployed from the microcatheter within the occluding thrombus, causing the thrombus embedded in the mesh of the stent. When it is pulled back, the stent mesh becomes deformed and grabs hold of the thrombus. The stent together with the thrombus is retrieved into the catheter and removed from the body. The objective of cerebral blood flow recanalization is thus achieved.


Figure 1 NeuroHawk™ thrombectomy device

Product Specifications/Use Guidelines

There are two models of thrombectomy device. Before determining that an thrombectomy device is to be used, it is necessary to estimate the diameter of the occluded vessel and the expected length of the occlusion site and to select the appropriate product model and compatible microcatheter according to the following table.

Table 1 Product models and specifications

Model	Nomi- nal diameter (mm) A	Stent working length (mm) B	Stent total length (mm) C	Length from fluoro-saver marker to distal end of stent (mm) D	Total length (mm) E	Delivery wire Max OD (inch) F	Minimum ID of compatible micro- catheter (inch)	Recom- ended vessel diameter (mm)
AIS4025	4.0	25	35	1500	2035	0.018	0.021	2.0 - 3.0
AIS6030	6.0	30	42	1500	2042	0.018	0.021	3.0 - 5.0

Intended Purpose

NeuroHawk™ thrombectomy device is intended to restore blood flow of brain vessels by removing thrombus for the treatment of acute ischemic stroke caused by large vessel occlusion.

Indication For Use

The product is indicated for the removal of thrombus in large vessels (including the M1 and M2 segments of the middle cerebral artery, internal carotid artery) in patients with ischemic stroke within 8 hours of symptom onset so as to restore blood flow.

Intended User

This product should be used by qualified professionals who have received necessary interventional training (especially for interventional radiology) to complete the endovascular procedures.

Target Population

NeuroHawk™ thrombectomy device can be used for patients who are diagnosed with acute ischemic stroke caused by large vessel occlusion (including the M1 and M2 segments of the middle cerebral artery, internal carotid artery), and need thrombectomy. No specificity in gender, race or nationality. Efficacy and safety in pediatrics, pregnant women, or children have not been evaluated.

Clinical Benefit

The benefits of the NeuroHawk™ Thrombectomy Device are to retrieve thrombus, achieve successful recanalization (mTICI, 2b–3) of previous occluded vessels, improve patients' functional independence (mRS 0-2), and reduce disability in the treatment of acute ischemic stroke due to large vessel occlusion.

Contraindications

- Patients with known allergies, resistance, or contraindication to antiplatelet agents, contrast agents, and/or anesthesia.
- Known allergy to Nitinol.
- Basic condition not suitable for endovascular intervention.
- Contraindicated for CT, CTA or MRI.
- Angiography suggests that the patient's vascular anatomy is not suitable for endovascular intervention. For example:

Severe intracranial vascular tortuosity;

Presence of intracranial vascular stenosis proximal to the obstructed vessel or implanted device that prevents the expansion or safe retrieval of the thrombectomy device;

Carotid artery dissection;

Other uncommon vascular morphology that may impede device use;
Inefficacy of drug therapy for intracranial arterial spasm.

Potential Complications

The use of this product may produce the following (but is not limited thereto):

- Headache
- Fever
- Nausea
- Vomiting
- Pain at puncture site
- Pseudoaneurysm
- Cardiac arrhythmia
- Hypervascularity/hypotension
- Myocardial infarction
- Infection
- Mental confusion
- Embolism
- Intracranial hemorrhage
- Hematoma
- Cerebral ischemia
- Hydrocephalus
- Neurological deficits
- Visual impairment/blindness
- Vascular dissection
- Perforation or rupture of a blood vessel
- Vascular spasm
- Vascular or tissue damage
- Intravascular thrombosis/vascular occlusion
- Adverse drug reactions, allergic reactions, including but not limited to: contrast agents, Nitinol and drugs
- Death

Precautions

- Use the product within the sterilization validity period.
- Store the product in dark, clean, and dry environment at room temperature with no corrosive gases.
- Do not expose the product to organic solvents (such as alcohol, etc.).
- Do not use if the inner packaging is found to be damaged.
- Before use, carefully inspect the sterile packaging as well as the system components to ensure that they have not been damaged during transportation. Do not use kinked or damaged components.
- Use the product in combination with a compatible microcatheter.
- The product should be operated under high quality x-ray observation. If resistance is encountered during manipulation, the cause of the resistance needs to be identified before proceeding with the surgical procedure.
- When using this product for thrombus removal, it is recommended, based on the product characteristics, that the number of removals for a single thrombectomy stent should be less than or equal to 5, and that the total number of removals in the same vessel should be less than or equal to 5 in order to reduce the risk of vascular injury.
- Select the appropriate thrombectomy device size so that its diameter after release is as close as possible to the diameter of the occluded vessel. If the wrong thrombectomy device size is selected, incomplete removal or injury to the vessel may occur.
- Care should be taken not to damage the distal end of the stent when introducing the thrombectomy device into the microcatheter from the introducer sheath.

Accessories Required for Product Use

- Microcatheter (inner diameter not less than 0.021 inch);

- Other accessories for performing a procedure to be selected based on the physician's experiences and preferences:
 - Guide catheter or long sheath (inner diameter not less than 0.053 inch);
 - Intermediate Catheters (5F or 6F);
 - Guidewire (outer diameter 0.035 inch, length 150 cm);
 - Micro-guidewire (outer diameter not greater than 0.014 inch, length not less than 185 cm);
 - At least two sets of normal saline or saline-heparin continuous flushing lines;
 - 60 ml syringe;
 - At least two units of rotating hemostasis valves;
 - At least two units of three-way stopcocks;
 - Infusion stand and femoral artery locking device.

Operation Procedures


Figure 2 Assembly diagram of thrombectomy device delivery

Preparation

- 1) Establish vascular access according to the standard operation method.
- 2) Determine the location of vascular occlusion and the diameter and length of the occluded vessel by angiographic fluoroscopy, and use antiplatelet and anticoagulant drugs according to medical standards.
- 3) Select the thrombectomy device product with the appropriate specifications according to the recommendations in Table 1. The effective length of the stent should completely cover the obstructed vessel thrombus.
- 4) Check the packaging of the corresponding product to ensure that the package is not opened or broken, and select the corresponding compatible microcatheter.

Warning: Improperly sized products may result in inability to reach the target lesion, low immediate recanalization rate, or vessel injury.

Placement of the guide catheter

- 1) Connect the RHV to the guide catheter, connect a 3-way stopcock to the side arm of the RHV, connect a continuous flushing device to the 3-way stopcock, and continuously flush the guide catheter by means of drip infusion.
- 2) Place the guide catheter into the target vessel using the standard method under the guidance of the guidewire, bring the distal end of the guide catheter close to the occlusion site, withdraw the guidewire from the guide

catheter, tighten the RHV, and maintain a continuous normal saline drip infusion.

- **Placement of microcatheter**

- 1) Connect the RHV to the microcatheter, connect the 3-way stopcock to the side wall of the RHV, connect the continuous flushing solution to the 3-way stopcock, and flush the microcatheter continuously by means of normal saline drip infusion.
- 2) Advance the microcatheter through the guide catheter to the target lesion using standard methods under the guidance of the micro-guidewire until the microcatheter passes beyond the thrombus (the distal end of the microcatheter should be placed beyond the distal end of the thrombus).
- 3) Withdraw the micro-guidewire from the microcatheter and tighten the RHV.

Warning: In case of resistance during the delivery of the product, the stent or the microcatheter should be withdrawn and checked for damage or kink.

- **Placement of thrombectomy device**

- 1) Open the package, fill the protective coiled tube with normal saline, remove the thrombectomy device from the coiled tube, push the stent out of the introducer sheath and carefully inspect it to ensure that there is no damage to the stent.
- 2) Resheath the stent into the introducer sheath.
- 3) Loosen the RHV on the microcatheter and insert the introducer sheath into the RHV but not fully into the microcatheter proximal hub, and tighten the RHV around the introducer sheath.

Warning: After the introducer sheath is inserted into the RHV of the microcatheter, ensure that the RHV on the microcatheter has been firmly tightened but not overtightened. Confirm that there are no bubbles in the system.

- 4) After waiting for flushing fluid to flow from the proximal end of the introducer sheath, loosen the RHV on the microcatheter again and advance the introducer sheath until it is fully seated within the microcatheter hub, and tighten the RHV but not so tight as to damage the stent during its insertion into the microcatheter;
- 5) Slowly advance the delivery wire so that the stent enters the microcatheter from the introducer sheath, and stop when the proximal end of delivery wire is less than 5 cm from the proximal end of introducer sheath;
- 6) Loosen the RHV on the microcatheter to remove the introducer sheath from the delivery wire. Put the introducer sheath within the sterile field, and do not discard it before the operation is finished;

Warning: To prevent gas from entering the blood, continue to deliver the thrombectomy device into the microcatheter only after there is fluid flow from the proximal end of the introducer sheath. Failure to withdraw the introducer sheath will result in interruption of flushing fluid injection and backflow of blood into the microcatheter.

- 7) Tighten the RHV on the microcatheter enough to advance the delivery wire but not much to compromise the continuous infusion. Continue to advance the delivery wire, turn on X-ray fluoroscopy after the proximal end of the RHV on the microcatheter is aligned with the fluoro-saver marks on the delivery wire.

Warning: The RHV must not be too loose or it will cause backflow of blood into the microcatheter.

- 8) Under fluoroscopic monitoring, slowly push the thrombectomy device within the microcatheter until the distal RO marker of the stent is aligned with the distal RO marker of the microcatheter (as shown in Figure 3, a); simultaneously adjust the microcatheter and thrombectomy device to ensure that the thrombus can be located in the middle-lower portion of the stent's working length after subsequent stent deployment.

Warning: If the thrombectomy device encounters abnormal resistance during pushing, stop pushing until the cause is found. When pushing the thrombectomy device inside the microcatheter, attention needs to be paid to the arrival position of the distal RO marker of the stent so that the stent is not unintentionally pushed out of the microcatheter, resulting in serious vascular injury.

- 9) Fix the delivery wire to prevent the thrombectomy device from migrating while slowly retracting the microcatheter in the proximal direction (as shown in Figure 3, b).
- 10) Continue to retract the microcatheter until the proximal markers of the stent emerge from the microcatheter, and stop retracting the microcatheter when the proximal markers of the stent are 5mm beyond the distal marker of the microcatheter (as shown in Figure 3, c-d).


Figure 3 Thrombectomy device deployment

Warning: After the thrombectomy device has been deployed, it cannot be pushed forward. Otherwise, the blood vessel would be seriously damaged. If the stent needs to be repositioned after deployment, the delivery wire should be fixed in position while pushing the microcatheter to retrieve the stent into the microcatheter. Then re-deploy it after the distal markers of the stent are flush with the distal marker of the microcatheter.

11) Keep the stent deployed for about 3 mins so that the thrombus becomes effectively embedded in the stent.

• **Retraction of the thrombectomy device**

- 1) Attach a 60 ml syringe to the side wall of the hemostasis valve of the guide catheter and continuously retract the syringe to aspirate the lumen of the guide catheter.
- 2) Adjust the position of the microcatheter until the proximal radiopaque coils of the stent are completely inside the microcatheter; tighten the RHV and secure the delivery wire and the microcatheter together, as shown in Figure 4.

Figure 4 Retrieval of the thrombectomy device

- 3) Maintain continuous aspiration of the guide catheter. Retract the thrombectomy device and the microcatheter all together simultaneously under negative pressure until the thrombectomy device and microcatheter are completely retrieved from the guide catheter.

Warning: If it is difficult to retract the microcatheter and thrombectomy device from the guide catheter, then the microcatheter, thrombectomy device, and guide catheter can be resheathed as a whole from the femoral sheath.

- 4) Observe the blood flow status of the occluded vessel by angiographic fluoroscopy. If the occluded vessel is effectively recanalized, the procedure is terminated by withdrawing the guide catheter and femoral sheath; if the occluded vessel is not effectively recanalized, a clinical assessment is made to determine whether to attempt another thrombectomy. If a second thrombectomy is required, it will be necessary to clean the stent, resheath it into the introducer sheath, and repeat the procedure in sequence from step 1) of the placement of the thrombectomy device until clinical efficacy is achieved.

Warning: When using this product for thrombus removal, it is recommended, based on the product characteristics, that the number of removals for a single thrombectomy stent should be less than or equal to 5, and that the total number of removals in the same vessel should be less than or equal to 5 in order to reduce the risk of vascular injury.

Disposal

After use, devices should be disposed of in accordance with medical waste disposal requirements.

Sterilization

The thrombectomy device has been sterilized using ethylene oxide and is non-pyrogenic. It cannot be re-sterilized for use.

Storage Conditions

Store in a dry, clean, well ventilated environment, at room temperature, with no corrosive gases.

Transportation Conditions

The product should be protected from stress, direct sunlight, rain or snows during transportation.

Shelf Life

When stored under the specified conditions, the thrombectomy device has a shelf life of 2 years. The specific production date is shown on the label.

Basic UDI-DI number

Neurohawk™ Thrombectomy Device Basic UDI-Device Identifier Number: 69586980T0048U4

Legal Statement

MicroPort NeuroTech (Shanghai) Co., Ltd. indicates definitely that the NeuroHawk™ Thrombectomy Device is intended for single use only. MicroPort NeuroTech (Shanghai) Co., Ltd. will not be liable for any incidental, special, or consequential loss, damage, or expense resulting, directly or indirectly, from improperly selecting size, inaccurate handling or any other man-caused accident. MicroPort NeuroTech (Shanghai) Co., Ltd. makes no warranty, representation or condition of any kind, whether expressed or implied (including any warranty of merchantability, suitability or fitness for a particular purpose) respecting the reuse of this stent. Moreover, MicroPort NeuroTech (Shanghai) Co., Ltd. assumes no responsibility or liability for incidental or consequential damages which may result from such re-use.

Meanwhile, MicroPort NeuroTech (Shanghai) Co., Ltd. assumes no responsibility or liability for any incidental or consequential damages which may result from the following circumstances:

- Any use, dismantling or disposal of the Product not in strict accordance with this Instructions for Use;
- Any intentional acts, gross negligence or non-professional acts of medical personnel;
- Damages caused by Force Majeure or accidents; and
- Other damages not caused due to the quality problems of the product

This document contains confidential information or proprietary technology of MicroPort NeuroTech (Shanghai) Co., Ltd.

After-sales Service

MicroPort NeuroTech (Shanghai) Co., Ltd. (hereinafter referred to as "MicroPort NeuroTech") takes "provide high quality, high-performance medical products" as the highest business purpose, and warrants that the entire product shall be free of defects in materials and workmanship upon receipt. Contact MicroPort NeuroTech directly with any questions related to the product.

Legal Manufacturer

Company Name: MicroPort NeuroTech (Shanghai) Co., Ltd.

Address: Building #16, 222 Guangan Road, Pudong New District, 201318 Shanghai, China

E-mail: Complaints-NeuroTech@microport.com

After-sales Service Company Information

Company Name: MicroPort NeuroTech (Shanghai) Co., Ltd.

English

Address: Building #16, 222 Guandong Road, Pudong New District, 201318 Shanghai, China
Tel: (86)(21)38954600
Fax: (86)(21)50801305
E-mail: Complaints-NeuroTech@microport.com

European Authorized Representative

MicroPort Medical B.V.
Paasheuvelweg 25, 1105BP Amsterdam, The Netherlands
Tel: +31 (0)20 545 0100 ext. 8
Fax: +31 (0)20 545 0109

If any device related serious incident occurs, please inform the Legal Manufacturer and/or European Authorized Representative and the competent authority of the Member State.

All rights reserved by MicroPort NeuroTech (Shanghai) Co., Ltd., and no reproduce without authorization of MicroPort.

“MicroPort” and logo“  ” are the registered trademarks of MicroPort® Sinica Co., Ltd., and are authorized to MicroPort NeuroTech (Shanghai) Co., Ltd., any other third-party without permission shall not use them.

Symbol Glossary

	Date of manufacture		Do not re-use
	Use-by date		Consult instructions for use
	Batch code		Contents 1 Unit
	Catalog number		Manufacturer
	Sterilized using ethylene oxide		Authorized representative in the European Community
	Do not use if package is damaged		Do not re-sterilize
	Keep away from sunlight		Unique device identifier
	Keep dry		CE Mark
	Single sterile barrier system		Medical device
	Caution		

More information on this device, including a Summary of Safety and Clinical Performance, is available at the websites below: <https://www.microportneurosci.com/>
<https://ec.europa.eu/tools/eudamed>

MicroPort NeuroTech (Shanghai) Co., Ltd.
Building #16, 222 GuangDan Rd.,
Pudong New District, 201318 Shanghai, China.
Tel: (86)(21)38954600
Fax:(86)(21)50801305
Email: Complaints-NeuroTech@microport.com
Website: <https://www.microportneurosci.com/>

MicroPort Medical B.V.
Paasheuvelweg 25
1105BP Amsterdam, The Netherlands
Tel: +31 (0)20 545 0100 ext. 8
Fax: +31 (0)20 545 0109

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Şirket Adı: MicroPort NeuroTech (Shanghai) Co., Ltd.
Tüm hakları MicroPort NeuroTech (Shanghai) Co., Ltd.'ye aittir, MicroPort'nin izni olmadan çoğaltılamaz.